


Mounting Instructions


SmartBucket seca 526 for seca 525

- Sicherheitshinweise • Safety precautions
- Consignes de sécurité • Avisos de seguridad


1.




2.




3.



- Lieferumfang • Scope of delivery
- Éléments livrés • Volumen de suministro





• Übersicht • Overview
• Aperçu • Vista general





seca 525

S/N > 1000000009050 ✓






S/N < 1000000009050 ✗




• Gerät ausschalten • Switching off the device
• Mettre l'appareil hors tension • Apagar el aparato




• Köcher demontieren • Removing the storage compartment
• Démonter le carquois • Desmontar el receptáculo




• SmartBucket montieren • Mounting the SmartBucket
• Monter le SmartBucket • Montar el SmartBucket




• USB-Verbindung herstellen • Establishing the USB connection
• Brancher le câble USB • Establecer la conexión USB



• Zubehör anbringen (variantenabhängig) • Attaching accessories (depending on version)
• Monter les accessoires (selon le modèle) • Colocar accesorio (dependiendo de la variante)




- Gerät einschalten • Switching on the device
- Mettre l'appareil sous tension • Encender el aparato




- SmartBucket/Temperaturmodul (variantenabhängig) aktivieren • Activating the SmartBucket/temperature module (depending on version) • Activer le SmartBucket/module de température (selon le modèle)
- Activar SmartBucket/módulo de temperatura (dependiendo de la variante)

1.



Login: admin
 Password: 1357

2.



Medical Measuring
Systems and Scales
since 1840

seca gmbh & co. kg
Hammer Steindamm 3-25
22089 Hamburg • Germany

Telephone +49 40 20 00 00 0

Fax +49 40 20 00 00 50

info@seca.com

All contact data under www.seca.com

seca
group

seca deutschland

seca france

seca united kingdom

seca north america

seca schweiz

seca zhong guo

seca nihon

seca mexico

seca austria

seca polska

seca middle east

seca brasil

seca suomi

seca américa latina

seca asia pacific